

**From:** [Rachel Freeman](#)  
**To:** [Sharon Reams](#); [Anja Doerge](#); [Lab Bend - All Lab Caregivers](#); [Lab Madras - All Lab Staff](#); [Lab Prineville – All Lab Caregivers](#); [Lab Redmond All Lab Caregivers](#)  
**Subject:** Re: GI PCR Testing Discontinuation for OP  
**Date:** Monday, January 13, 2025 11:24:28 AM  
**Attachments:** [image.png](#)

---

Hi All -

I've got an IT ticket in to update the containers for the **NOROVIRUS (NORWALK) PCR [LAB495]** test, Mayo test code LNORO. The preferred container for this test is the **orange top Para-Pak C&S vial**. It's currently set to default of a standard Cary Blair, which is also accepted as an alternate container.

**C and S Vial (T058)**



However, Mayo does not accept the yellow-top MCC Modified Cary Blair for this test and will result in cancellation.



Thank you,

**Rachel Freeman** | Referral Specialist II | Work: 541-706-3722 | Fax: 541-706-6326 | [freeman.rachel@mayo.edu](mailto:freeman.rachel@mayo.edu)

**Mayo Clinic Laboratories** | 3050 Superior Drive NW | Rochester, MN 55901 | Toll Free: 800-533-1710 | [www.mayocliniclabs.com](http://www.mayocliniclabs.com)

---

**From:** Sharon Reams <[sreams@stcharleshealthcare.org](mailto:sreams@stcharleshealthcare.org)>  
**Sent:** Thursday, December 12, 2024 9:27 AM